BIO SYNEX

Autotest COVID-19

Test destiné à l'autodiagnostic.

INTRODUCTION

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Test for self-diagnosis.

INTRODUCTION

COVID-19 is an acute respiratory illness caused by infection with SARS-CoV-2. SARS-CoV-2 belongs to the same family of viruses that cause Severe Acute Respiratory Syndrome (SARS). The most common symptoms of COVID-19 are fever, fatigue and a dry cough. Some patients may experience pain, nasal congestion, runny nose, sore throat, diarrhea, and loss of taste and smell. These symptoms are usually mild and begin gradually. Some people are infected but do not develop any symptoms. Most people (about 80 %) recover from the disease without the need for hospital care.
The disease can be spread from person to person through small droplets emitted from the nose or mouth when a person infected with COVID-19 coughs or exhales. Transmission through contact of contaminated hands with the face is common. Hands can become contaminated after touching surfaces that carry the virus. The estimated incubation period for COVID-19 ranges from 1 to 14 days.

INTENDED USE

COVID-19 self-test is a rapid immunochromatographic self-test for the qualitative detection of IgG and IgM antibodies against SARS-CoV-2 from a drop of capillary blood. The test comprises SARS-CoV-2 antigens (called RBD - Receptor Binding Domain) coupled with colloidal gold particles that are capable of detecting IgG and IgM antibodies directed against SARS-CoV-2. This self-test determines whether you have been in contact with the virus responsible for COVID-19. The COVID-19 self-test, based on antibody detection, can identify individuals who are unaware that they have been infected, either because they have never developed symptoms or because they have symptoms that have never been properly diagnosed. This means that self-testing can identify silent infections, as well as people who have been sick but have recovered. This test is for use only in adults or under adult supervision if performed on a child.

The COVID-19 self-test will give a positive result in the majority of cases after 13 days of the onset of symptoms. A negative result with COVID-19 Test in a symptomatic individual does not exclude COVID-19 infection. If symptoms suggestive of COVID-19 are present, a negative result should be followed up with a second test 11 days after the first test.
This self-test does not detect the presence of the virus responsible for COVID-19, but it does detect the response of a patient's immune system against the COVID-19 virus.

KIT CONTENTS

*3 pipettes are provided in case of loss or breakage. Ref 859082 contains one bag for one self-test.

TESTING GUIDELINES

- Wash your hands with warm water beforehand and dry them properly. Warm water improves circulation at the fingertips.
- You can let your hand hang for about a minute along your body so that more blood can flow into it. It is recommended that you prick your left hand for a right-handed person and your right hand for a left-handed person.
- It is best to prick the fingertip on the side: this is where there is the most blood and less sensitivity. Avoid areas where the skin is thicker and therefore more difficult to prick.
- It may happen that a drop does not form right away. In this case, you should wait a while and massage your finger lightly, from the base to the puncture site.
- Do not hesitate to form large drop of blood at the tip of your finger. This way, the pipette will fill up with blood easily and quickly.

LIMITATION OF THE TEST

- This test is not a PCR screening test to detect the presence of the virus, but it does detect a reaction of the patient's immune system against COVID-19.
- This test can be used to identify silent infections with no symptoms, past or present, and people who have been sick but have recovered.
- This test does not determine whether the patient is contagious, it only indicates whether the person has been in contact with the virus at some point.
- A false negative result may occur in some patients with early symptoms. Repeat the test on subsequent days to detect the development of antibodies. At the onset of infection, anti-SARS-COV-2 IgM concentrations may be below detectable levels.
- The test will only indicate the presence of anti-SARS-COV-2 antibodies in the sample and should not be used as the sole diagnostic criteria.
- This test is not a substitute for a medical consultation or for the result of a biological analysis performed in a medical laboratory.
- Incorrect results may be obtained:
 - If the test is not used according to the instructions in the package insert,
 - If the foil pouch is damaged or if the test is not performed immediately after opening the foil pouch,
 - If the storage conditions are not met or if the test is performed after the expiry date on the foil pouch.
- A positive result must be confirmed by laboratory analysis. Consult your doctor.
- There is a possibility of obtaining false positive results, especially in some cases of patients with rheumatoid arthritis.
- Cases of cross-reactions have been observed with SARS-CoV-1 and other seasonal coronaviruses.
- If you have a positive result showing the presence of antibodies, it does not mean that you are immune to COVID-19.
- The response of the immune system is unique to each individual. Antibodies may disappear in some individuals as soon as a few weeks after infection.

WARNINGS

- Make sure that the blood has been deposited before applying the drops of diluent.
- Keep the self-test and its components out of the reach of children; the diluent can be dangerous if swallowed and the lancet a risk of injury.
- The used test cassette as well as all components of the test can be disposed of with household waste in a well-sealed bag.
- All components included in this kit are for this test only. Do not reuse this test and/or its components.
- This test should be performed at room temperature (15°C – 30°C).
- This test is for single use only.
- This test should be used immediatly after collection of the blood drop from the fingertip.

STORAGE

- This test should not be frozen and should be stored in a dry place between 2°C and 30°C.
- The test cassette should remain in the sealed pouch until use.
- Do not use the self-test beyond the expiration date indicated on the foil pouch and on the outside of the carton.

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The result of the test is **NEGATIVE** if a colored line appears only in the control zone (C) and no colored line appears in the IgG and IgM zones (regardless of the intensity). A negative result indicates the absence of antibodies against COVID-19. Your body has not been in contact with COVID-19. A negative result does not rule out recent SARS-CoV-2 infection. If you think you have been in contact with the virus (with an infected person) in the last few days before the test, we invite you to do a second self-test 10 days after your last exposure. Continue to protect yourself and keep a social distance.

EN The result of the test is **POSITIVE** if 2 or 3 colored lines appear on the membrane. One colored line appears in the control zone (C) and one in the IgG and/or IgM zones, regardless of the intensity of the lines. A positive result (cases 1, 2 or 3) shows that your body has been in contact with the COVID-19 virus and that your immune system has produced IgM and/or IgG antibodies.

A positive result with the Biosynex COVID-19 self-test does not tell you whether you are still a carrier of the virus and therefore a transmitter.

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