



Lab Quality Results at the Point of Care¹

The A1CNow⁺ test system provides healthcare professionals with a fast, easy, and accurate method of obtaining HbA1c results anywhere, including the exam room. In 5 minutes, clinicians have information for real-time consultation with a patient.




A1CNow⁺ Test System Key Features

- Lab-quality accuracy
- Portable
- Easy to use
- Fingerstick (5 μ L) blood sample
- NGSP-certified
- IFCC-traceable
- CLIA-waived

The A1CNow⁺ test system has demonstrated accuracy equivalent to clinical laboratory analyzers!¹

Comparison of A1CNow⁺ Relative to the Clinical Laboratory¹



Accuracy

Average Bias = -0.03 %A1C
Bland-Altman (95% confidence) =
-0.6 to +0.5 %A1C

Lab Analyzers Used in Evaluations:

Tosoh G8,
Roche Cobas® Integra 400 plus,
Roche Cobas c513,
Roche Cobas Integra 800,
Beckman Coulter AU5400,
Siemens Dimension®

Precision

(14 samples / 10 replicates each)
Average %CV = 3.34%

NGSP Certification Ensures Accuracy

The A1CNow⁺ analyzer was designed to be as accurate as sending blood to a lab, allowing medical professionals to have informed conversations about how patients are managing their diabetes in minutes, not days. The A1CNow⁺ test system is annually certified by the National Glycohemoglobin Standardization Program (NGSP).

Available
in a 10-count
or 20-count
test kit



Sold
Separately



1. Summary of 2017 Evaluations Comparing the A1CNow⁺ System to Clinical Lab Analyzers by PTS Diagnostics. Data on file.



Website: ptsdiagnostics.com/a1cnow-accuracy/
Office: +1 (317) 870-5610 FAX: +1 (317) 870-5608
Toll-Free in the US: 1 (877) 870-5610