


ez-Track¹



ez-Track¹

The information contained herein is subject to change without prior notice.
Theradiag shall not be liable for any technical and editorial errors or omissions contained herein.

Additionally, any malfunction in elements including the optical sensor assembly, the system board and the drive train should be consulted with and serviced by Theradiag.
Any unauthorized operations on elements not suggested in this document will void the existing warranty.

Any serious incident that has occurred in relation to this product shall be informed to Boditech Med Inc., and your competent authority.

For technical assistance, contact:

Theradiag's Technical Services

Tel: +(33)-1-64-62-10-12 / Fax : +(33)-1-64-62-09-66

E-mail : support@theradiag.com



Boditech Med Inc.

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,
Gang-won-do 24398, Republic of Korea
Tel: +(82)-33-243-1400 / Fax: +(82)-33-243-9373
www.boditech.co.kr



Theradiag

14 Rue Ambroise Croizat 77183 Croissy Beaubourg France
Tel: +(33)-1-64-62-10-12 / Fax : +(33)-1-64-62-09-66
www.theradiag.com



Obelis s.a.

Bd. Général Wahis 53, 1030 Brussels, Belgium
Tel: +(32)-2-732-59-54 / Fax: +(32)-2-732-60-03
E-Mail : mail@obelis.net



The actual contents of the display could differ, depending on the specification and customer requirements.











Table of Contents














1. General Warnings and Important Instructions	3
2. Principle of Operation	5
3. Commercial Package Contents	9
4. Specifications	11
5. Operational and Functional Elements	12
6. Installation	15
7. Running Tests	19
8. Supportive Procedures	30
9. Configuring	50
10. Maintenance, Servicing and Disposal	72
11. Trouble shooting	79
12. Warranty	82

1. General Warnings and Important Instructions

Graphic symbols

The following list explains various symbols, which the user will find in this operation manual as well as in package inserts and on labels of ez-Track¹ products.

Symbol	Description
	Manufacturer
	Date of manufacture
	Authorized representative in the European Community/European Union
	CE conformity marking
	For <i>in vitro</i> diagnostic use
	Catalogue number
	Serial number
	<p>Biohazard</p> <p>To warn of a biological hazard. Taking care to avoid exposure to a biological hazard.</p>
	<p>Biological risks</p> <p>To indicate a reference to substances that may be hazardous to men, animals, plants, or the environment based on biological activity (for example, holding a virus).</p>
	<p>Caution</p> <p>This symbol is essentially a safety symbol and should be used to highlight the fact that there are specific warnings or precautions associated with the device, which are not otherwise found on the label. The symbol "Caution" is still sometimes used to have the meaning of "Attention, see instruction for use".</p>

	Indicates a hazardous situation, which if not avoided, could result in injury to the operator or a bystander (e.g. electrical shock or UV exposure).
	LASER Radiation Do not stare into beam
	Consult instructions for use
	Separate collection for electrical and electronic equipment
	USB memory
	Temperature limitation
	Humidity limitation
	Direct current
	Alternating current
	On (Power)
	Off (Power)
	Importer
	Distributor



Danger!

Improper use or careless handling of ez-Track¹ may cause electrocution, burns, fire and other hazards. To reduce the risk of electric shock:

- Use only the adaptor provided with ez-Track¹.
- Never try to disassemble ez-Track¹ as it does not have any operator-serviceable component.
- Do not allow ez-Track¹ to come in contact with water or any liquid in any manner.
- Only occasional cleaning of cartridge holder and exterior by wiping with a soft cloth is recommended. Unplug ez-Track¹ from external power source before cleaning.
- Disconnect the adaptor during extended period of non-use and/or storage.



Laser Radiation

ez-Track¹ has an internal laser module. It can cause some risks due to laser radiation. To protect these risks, please keep the followings.

- Do not alter or disassemble the ez-Track¹. Unauthorized alternation and disassembly can cause exposure to laser radiation. When operating the ez-Track¹, please keep and follow the instructions in this operation manual.
- There is no need to use specific eyes protective gears because it is impossible to see laser beam during the normal use.
- Do not open the holder cover during a test. If you look internal laser beam by using a microscope or telescope, it may cause some risks.
- The exteriors of the ez-Track¹ protect users from laser beam. Do not open the exteriors.
- It is impossible to adjust the laser power except the manufacturer.



Biological risk

Patient samples and used ez-Tracker cartridges must be treated as being potentially infectious material. To reduce the risk of any biohazard:

- Used cartridges should be disposed of in accordance with local regulations as applicable to the disposal of bio-hazardous materials.
- At the end of its useful life, ez-Track¹ could be sent back to Theradiag. for recycling or proper disposal. Alternatively, it should be disposed of in accordance with relevant local regulations.
- To reduce the risk of any biohazard or chemicals, wear protective gloves. Do not touch or use test cartridges and holder cover with bare hands.

Always dispose of biological waste in accordance with relevant regulations.



Caution

ez-Track¹ is meant for countertop operation under laboratory conditions.

To reduce the risk of damage to ez-Track¹:

- Handle ez-Track¹ gently with care to prevent damages to its internal components.
- Keep ez-Track¹ in a dust-free place.
- Do not expose ez-Track¹ to an environment having high humidity, direct sunlight and/or extreme temperature.
- Do not place anything on top of ez-Track¹.
- Use only the accessories provided with ez-Track¹ and cartridge(s) and ID chip(s) supplied by Theradiag.
- Unless otherwise instructed specifically, do not insert anything into any opening, socket or drawer provided in ez-Track¹.
- During transportation or storage, keep ez-Track¹ in dry conditions having -20~50 °C (-4~122 °F) temperature in an upright position and protect it from mechanical shocks.
- Before moving to another location, remove tips from the tip box.

Failure to follow these precautions will invalidate the warranty.

To reduce the risk of obtaining misleading or erroneous results with ez-Tracker tests:

- Ensure that ez-Track¹ is in proper working condition before performing the test.
- ez-Tracker test setup should be arranged in an operational environment having temperature 15~35 °C (59~95 °F) and humidity 10 ~ 70 %.
- Use ez-Tracker cartridges that have been stored in recommended storage conditions only.
- Expired cartridges must not be used.
- Before performing any ez-Tracker test, ensure that name of the test as well as lot number printed on the 'ez-Tracker ID chip' exactly match with those of the 'ez-Tracker cartridge'.
- Do not remove the cartridge from its blister packing until testing and allow it to stand at room temperature for at least 30 minutes.
- Before performing any ez-Tracker test, read the entire package insert of the test thoroughly and follow the instructions carefully.

Failure to follow these precautions may yield misleading or erroneous results with the ez-Track¹ test.

2. Principle of Operation

Intended Use

ez-Track¹ is an analyzer intended for use in conjunction with fluorescence immunoassay (FIA) kits for quantitative, semi-quantitative and qualitative measurements of various analytes.

For *in vitro* diagnostic use only.

Intended User

ez-Track¹ should be operated by healthcare professionals or trained medical technologists (professional use only).

Principle

ez-Track¹ is a fluorescence and europium nanoparticle scanning instrument to be used in conjunction with various ez-Tracker tests which are based on antigen-antibody reaction and fluorescence technology (Fluorescence and europium nanoparticle immuno-assay).

ez-Track¹ uses a semiconductor diode laser as the excitation light source for illuminating the sample-loaded ez-Tracker cartridge(s) inserted in its cartridge slot(s); thereby triggering fluorescence from the fluorochrome-labeled complexes of analyte accumulated at the test line on the cartridge membrane.

Fluorescent light is collected together with the scattered laser light. Pure fluorescence is filtered from the mixture of scattered and fluorescent light. Intensity of fluorescence is scanned and converted into an electric signal which is correlated with the concentration of target analytes in the sample.

The on-board microprocessor computes concentration of the analyte in test sample based on a pre-programmed calibration derived from the 'ez-Tracker ID Chip' inserted in the ID chip port.

Computed and converted test result is displayed on the display screen of the ez-Track¹ in a quantitative manner.

ez-Tracker tests are *in vitro* diagnostic tests intended for screening, monitoring and/or routine physical examination.








3. Commercial Package Contents

ez-Track¹ and its accessories are provided in a single box. For proper operation of the system, user must be familiar with the technical term and use of each component.

User must ensure that following items are included the commercial package. If any item(s) were to be missing, please contact your sales distributor or Theradiag’s support clients at +(33)-1-64-62-10-12 in France.




■ **Standard items**

1. 1 ez-Track¹
2. 1 Operation Manual
3. 1 Set of ez-Track¹ System Check Cartridge and System Check ID Chip
4. 1 Set of AC adaptor and power cable
5. 1 Thermal printer paper roll
6. 1 ez-Track¹ OS SD Card (Inserted in the instrument)

1	2	3a
		
3b	4	5
		
6		
		

■ **Optional items**

1. Barcode Scanner
2. USB Keyboard (86 keys)
3. USB OTG Cable

1	2	3
		

4. Specifications

Physical Features

<input type="checkbox"/> Dimensions	320.6 mm (L) x 180.6 mm (W) x 206.4 mm (H)
<input type="checkbox"/> Weight	4.0 Kg
<input type="checkbox"/> Power supply	DC 12V / 5A AC/DC Adaptor Input : 100-240 V AC, 50~60 Hz, 1.5 A Output : DC 12V / 5A
<input type="checkbox"/> Data output	LCD / Internal printer / USB OTG port

Environmental Requirements

<input type="checkbox"/> Operation Temperature	15 ~ 35 °C (59 ~ 95 °F)
<input type="checkbox"/> Storage Temperature	-20 ~ 50 °C (-4 ~ 122 °F)
<input type="checkbox"/> Humidity	10 ~ 70 %
<input type="checkbox"/> Location	Dry, clean, flat and horizontal surface in an indoor setting away from direct sunlight, mechanical vibration and magnetic field

Optical Components


<input type="checkbox"/> Light source	Laser Diode
<input type="checkbox"/> Detector	Silicon Photo Diode

Other Components

<input type="checkbox"/> Interface	USB OTG	1 Port (only for manufacturer)
	USB	4 Ports
	LAN	1 Port (only for manufacturer)
	SD Card Slot	1 Port (only for manufacturer)
<input type="checkbox"/> Display	LCD (800x480 Graphic LCD, Touch-screen module)	

ez-Track¹ complies with the EMC requirements according to EN 61326-2-6.
ez-Track¹ complies with the safety requirements according to IEC 61010-1 and IEC 61010-2-101.

5. Operational and Functional Elements



No.	Description
① LCD display	Acts as the platform for inputting all information regarding the test. Displays test result and message alerts.
② Power Button	For switching ez-Tracker ¹ 'On' or entering it in 'Sleep Mode'
③ ID Chip port	For inserting ID Chip
④ Holder	Accommodates cartridge holders meant for inserting ez-Tracker test cartridges
⑤ Tip box	Collects used pipette tips
⑥ Internal Printer	Offers print-out of test result

⑦ Power port	For connecting the AC adaptor
⑧ USB OTG port	For connecting PC (only for manufacturer)
⑨ USB port	For connecting Barcode scanner / Keyboard
⑩ Main Power switch	For switching ez-Track ¹ on/off
⑪ LAN port	For connecting to the internet (only for manufacturer)
⑫ SD Card port	For software upgrade

Function key




Power Button

If you press this button momentarily, ez-Track¹ will enter/exit 'Sleep mode' and the LCD screen will be turned off/on.


If you press this button for more than 4 seconds, ez-Track¹ will shut down.

5.1 Software menu tree

Menu tree of ez-Track¹ is as follow.



6. Installation

 Important	Before using ez-Track ¹ , read the warnings and precautions described previously in 'Clause 1: General Warnings and Important Instructions' of this manual.
---	--


Power supply

ez-Track¹ will work on household power line i.e. 100-240 V AC, 50~60 Hz. If you are uncertain of your power line specifications, consult your utility manager.

Operational environment


ez-Track¹ should be placed on a dry, clean, flat, and horizontal surface in an indoor setting away from direct sunlight, mechanical vibration and any source of electromagnetic radiation (e.g. power transformers). Room temperature of the operational chamber must be maintained in the range of 15~35 °C (59~95 °F).

Installation procedure

	Remove the contents from the package and inspect the contents for any physical damage and/or missing item(s).
---	---

ez-Track¹ is provided with a foam for protecting the inside of ez-Track¹. When one turns on the ez-Track¹, the holder will come out with a foam.



	Please remove this foam before using the ez-Track ¹ as shown above.
---	--

1. Connect the AC adaptor by inserting its connector into the power port.



2. Connect the power cord to external power outlet and turn on the main power switch.



3. Press the power button provided on the front of ez-Track¹.



4. Upon pressing the power button, ez-Track¹ will display Theradiag's logo and it executes a series of system initialization and self-checking routines.




During system initialization, ez-Track¹ will display each test items, and the test results will be displayed. If there is any problem, an error message is displayed. If a user taps the "Confirm" button, he/she can use ez-Track¹ except for test functions. If you want to perform tests, please turn off ez-Track¹ using the main power switch (provided on the rear side of ez-Track1) and restart ez-Track¹. If the same problems are present, please contact the local distributor or Technical Service of Theradiag.



5. After finishing the initialization, ez-Track¹ prints out "System OK!" and displays the popup window for selecting registered users.
6. When you select a user, ez-Track¹ displays "Log in" window. If you want to change the current user to another user, please tap the drop-down button and change the user.

7. Running Tests

	<p>Before performing any ez-Tracker Test, please read the respective package insert carefully for other details of the test procedure.</p>
	<p>Before performing any ez-Tracker Test, ensure that name of the test as well as the lot number printed on the 'ez-Tracker ID Chip' exactly match with those on the 'ez-Tracker Cartridge.' Otherwise ez-Track¹ will display an error message "Insert ID chip".</p>
	<p>ez-Tracker Cartridges and ez-Tracker Controls must have been stored under recommended storage conditions and must not have passed the respective expiration date otherwise, the test may yield misleading/ erroneous results.</p>
	<p>If scanning of cartridge is terminated abruptly or any problem occurs, a relevant error code or a message alert will be displayed. Refer to 'Clause 11: Troubleshooting' for probable cause(s) and remedial/corrective measure(s) for the displayed error code or message alert.</p>

7.1 Test Setup


- For performing an ez-Tracker test, arrange the following items in addition to ez-Track¹ in a test-ready mode:
- ez-Tracker XXX ID Chip(s)
- ez-Tracker XXX Cartridge(s)
- Patient Sample(s) or ez-Tracker XXX Control(s) or Calibrator(s)
(Note: XXX represents name of the ez-Tracker Test being performed.)
- Ensure that the name of the test as well as the lot number printed on the 'ez-Tracker ID Chip' with those printed on the 'ez-Tracker Cartridge' before removing the cartridge from its blister packing.
- Store patient sample(s) or control(s) and the cartridge(s) at room temperature for about 30 minutes before testing.

7.2 Logging in and Inputting required information

- As described in 'Clause 6: Installation', test-ready ez-Track¹ displays the "Log in" window.
- Input user name and password. The default 'user name' and 'password' are "admin" and "0000," respectively. You can change the password for 'admin' and add other users. Please refer to 'Clause 9.2: User management; for further details. If log in is successful, main menu will be displayed.

7.3 Main menu description

7.3.1 Details of the information/main menu displayed on the screen are as follow.



The screenshot shows the main menu of the Theradiag ez-Track device. At the top, there are three tabs: 'Analyses', 'History', and 'Cal/QC'. To the right of these tabs is a hamburger menu icon. Below the tabs is a 'Patient ID' input field. Underneath are fields for 'Age (Optional)' and 'Gender (Optional)' with radio buttons for 'Male' and 'Female'. There are two sections for 'Tip type' and 'Sample type'. The 'Tip type' section has 'General' and 'C-tip' options. The 'Sample type' section has 'Whole blood', 'Serum/plasma', and 'Others' options. A large box labeled 'ITEM' contains 'Reagent Lot' and 'Time remaining'. Below this is a status indicator 'Idle'. At the bottom, there are four large buttons: a green 'Waste tip' button (0%), a grey 'Holder In/Out' button, a red 'Stop' button, and a blue 'Start' button. The bottom status bar shows the user 'admin', a settings gear icon, a LAN connection icon, and the date/time 'Aug-15-2016 / 22:04' and version 'AA1.VW.01.01.23'.

Extended menu:

- System
- User management
- ID chip Management
- Temperature

QC

Check, search and print previous test result(s).

▪ Test information
 ▪ Remaining reaction time

Holder In/Out

Waste tip chamber status indicator. It has "Reset" function.

Log in / out

'User Name' currently logged in.


Setting for QC

LAN connection status indicator
 White : LAN is connected.
 Grey : LAN is disconnected.

Stop the current test

Start the test

▪ Current time
 ▪ S/W version number




7.3.2 When you start the test after inserting ID Chip, information of that ID Chip will be stored/registered in ez-Track¹ automatically. If you want to test the same lot of ez-Tracker test, ID Chip is not needed.

To register ID Chip manually, refer to 'Clause 9.3 ID Chip Management'.

7.4 Inputting test information

7.4.1

- **Inputting required test information'**
 - Select "Tip type" and "Sample type".
 - If you select "General" as tip type, you should select "Sample type" among Whole blood, Serum/Plasma and Others.
 - If you select "C-tip" as tip type, the sample type is selected as "Whole blood" automatically.



7.4.2

■ **Optional information**

- It is possible to perform a test without inputting the following 3 information.

■ **Inputting Patient ID**


- Input a patient ID using a keyboard which is displayed on the LCD screen. When you tap the 'Done' button, the keyboard will disappear. If USB keyboard is connected, you can use it to input a patient ID.

■ **Inputting Age**



- Input patient's age using a numeric keypad which is displayed on the LCD screen. When you tap the 'Done' button, the keypad will disappear. If USB keyboard is connected, you can use it to input patient's age.

■ **Gender**



- Tap the blank box for 'Male' or 'Female'. If you want to cancel the selection, tap the selected box again.



7.5 Inserting Sample-loaded Cartridge and Obtaining Test Result


<p>7.5.1</p>	<ul style="list-style-type: none">■ Considering that 'ez-Tracker XXX ID Chip' has been inserted in the ID chip port, proceed as follows for performing the ez-Tracker XXX test on clinical sample:<ul style="list-style-type: none">○ Load recommended volume of the test sample into the 'Sample Well' of 'ez-Tracker XXX Cartridge'. <p>(⚠ Refer to the package insert of the ez-Tracker test being performed for details on sample processing, volume to be loaded, etc.)</p>	
<p>7.5.2</p>	<ul style="list-style-type: none">■ After inserting the test cartridge into the specific cartridge holder, place a 'Pipette Tip' (with its broad end facing upwards) into the 'Tip Slot' of the test cartridge.■ While inserting the cartridge, sample well should be oriented outward from ez-Track¹. <p>(⚠ For inserting the cartridge, push it gently with your thumb or index finger; all the way into the cartridge holder until it comes to a stop. Avoid applying excessive force; otherwise it may result in mechanical failure of the system.)</p>	





<p>7.5.3</p>	<ul style="list-style-type: none"> ■ To proceed with the test, tap the “Start” button. ■ The new window pops up. Please check the inputted test information (Tip type / Sample type / Patient ID / Age / Gender). ■ To proceed with the test, tap the “Confirm” button. ■ If you want to cancel the test, please tap the “Cancel” button. 	
<p>7.5.4</p>	<ul style="list-style-type: none"> ■ ez-Track¹ reads the barcode of inserted ez-Tracker cartridge and check the test item and the lot number automatically. ■ ez-Track¹ checks whether a pipette tip is inserted or not automatically. ■ The name of test and lot number will be presented in the LCD display(refer to ①). ('PCT' is the name of the test and 'PCLXB03F' is the lot number.) ■ While the test is in process, <u>“In progress”</u> is displayed(refer to ③). 	

	<ul style="list-style-type: none"> The remaining time is displayed in minutes and seconds (mm:ss) on the right of the test name(refer to ②). 	
<p>7.5.5</p>	<ul style="list-style-type: none"> It is impossible to start a test if a pipette tip is not inserted or a test cartridge is not inserted completely. In this case, ez-Track¹ displays a “Warning” window and the test does not proceed. You need to solve the problems(insert a pipette tip, or a test cartridge) and tap ‘OK’ button to start a test. For proper operation, tap the ‘Start’ button, then in a new popup display, tap ‘Confirm’ to proceed with further steps. 	
<p>7.5.6</p>	<ul style="list-style-type: none"> If a test cartridge is not inserted, ez-Track¹ displays ‘Barcode error’. If ez-Track¹ displays ‘Barcode error’ although a test cartridge is already inserted, please remove and insert the test cartridge again. Then proceed with the test. If ez-Track¹ displays ‘Barcode error’ continuously, please contact Client Support center of Theradiag or regional distributor. You need to insert a test cartridge and resume the test procedure according to clause 7.5. 	





7.5.7




- When the test is completed, ez-Track¹ will display a new pop-up window. You can find the test information (mandatory and optional information), cartridge information and the test result. In addition, ez-Track¹ will print out the test result using the internal thermal printer.
- ez-Track¹ can save 5,000 test results in the internal memory. (To check the saved test results, please refer to clause 8.2.)



	<p>Before performing any ez-Tracker Test, ensure that the name of the test as well as the lot number mentioned on the 'ID Chip' must match with that of the 'ez-Tracker Cartridge'.</p>
	<p>'ez-Tracker Cartridge' having a lot number mismatching with that of the 'ID Chip' will lead to display an error alert, 'Insert ID chip'.</p>
	<p>The used cartridge should be disposed of in accordance with local regulations regarding the disposal of bio-hazardous materials.</p>
	<p>If you press the power button during the test, the display will be turned off. However, the test run will proceed normally without affecting the test result.</p>

7.6 Finishing Working With ez-Track¹

<p>7.6.1</p>	<ul style="list-style-type: none"> Before shutting down ez-Track¹, ensure that ID chip(s) and cartridge(s) have been removed from ID chip port(s) and cartridge holder(s), respectively. 	
<p>7.6.2</p>	<ul style="list-style-type: none"> Press the 'Power Button' for more than 4 seconds. 	
<p>7.6.3</p>	<ul style="list-style-type: none"> 'Power off' menu will be displayed. Tap 'OK' to shut down ez-Track¹. 	
<p>7.6.4</p>	<ul style="list-style-type: none"> 'Power off' will be executed. The cartridge holder will retract automatically. ez-Track¹ will shut down. 	

7.6.5	<ul style="list-style-type: none">■ Turn off the main power switch provided on the rear side of ez-Track¹.	
7.6.6	<ul style="list-style-type: none">■ If you do not intend to use the system for a long period (e.g. more than 2 days), unplug ez-Track¹ power cord from the external power outlet.	
7.6.7	<ul style="list-style-type: none">■ Remove the power cable from the 'Power Port'.	

8. Supportive Procedures

For ensuring efficient and reliable performance as well as proper recording and reporting of test results, the overall operation of ez-Track¹ essentially involves following supportive procedures in addition to running ez-Tracker tests.

8.1 System check, Calibration and QC

- System check
 - To perform a system check, a user needs to use a specific system check cartridge and a system check ID chip.
 - ez-Track¹ measures the intensity of fluorescence and compares that result with pre-set value to inspect if the measured signal agrees with pre-set values.
 - ez-Track¹ repeats the measurement several times to ensure stable performance of a fluorescence system when checking the system.



- The system check is needed as below:
 - When ez-Track¹ is installed initially.
 - When you move ez-Track¹ to another place and install it again.
 - When some errors are generated.
 - When a test result is considered abnormal or differs from the reference value.





- Calibration
 - The master calibration curve for cartridges is stored in the ID chip included in the cartridge box and no separate calibration process is required for the new reagent.
 - Calibration curves stored in the ID chip may not be valid for all test environments.
 - In this case, it is required to carry out the calibration process.
 - A dedicated calibrator and a test cartridge are required for calibration.

- The calibration is needed as below:
 - When the QC test results are out of range.
 - When the test results are suspected to be abnormal.
 - ※ ez-Tracker test system is designed for point of care, emergency, and small-volume testing. Therefore, a regular calibration is not required/recommended.

- QC (Quality control)
 - To perform QC, a user needs to use QC reagent which is supplied by the manufacturer on demand.


- The QC is needed as below:
 - When a test lot is changed.
 - When a test result is considered abnormal or differs from the reference value.
 - When a test result is not consistent with symptoms.
 - For teaching or learning the system.

<p>8.1.1</p>	<ul style="list-style-type: none"> ■ To do calibration and QC, you need to activate this function. ■ Select the ① 'Cal/QC' tap. ■ To activate 'Cal/QC' function, please tap the setting icon ②. 	
<p>8.1.2</p>	<ul style="list-style-type: none"> ■ Tap the 'Enable Calibration' menu. The check box will be toggled whenever you tap this menu. ■ To activate the calibration and QC function, the check box should be checked with <input checked="" type="checkbox"/>. 	

<p>8.1.7</p>	<ul style="list-style-type: none"> ■ 'Add control' window will pop-up. ■ Input a lot number, item name, mean value, low value, high value and expiration date using key pad or a barcode reader. <p><i>Note: The expiration date is set as YY-MM. The actual date will be set as the last day of the selected month. For example, if you set expiration date as 2015. 11, the detailed expiration date is 2015. 11. 30.</i></p> <ul style="list-style-type: none"> ■ To apply this information, tap the 'OK' button. 	
<p>8.1.8</p>	<ul style="list-style-type: none"> ■ If you want to delete a registered control, select a control that you want to delete. ■ Tap the  icon to delete a registered control. ■ Tap the 'OK' button for confirmation. 	
<p>8.1.9</p>	<ul style="list-style-type: none"> ■ To do calibration and QC, please prepare a test cartridge that you want to test. ■ Dispense a calibrator and control into the sample well of the cartridge. ■ Insert the cartridge into the cartridge holder and insert tips. 	



8.1.10

- Tap the 'Cal 1' button to perform Cal/QC.




8.1.11

- When you tap 'Cal 1' button, 'Cal-1 information' window pop-up.
- You can select item and lot number that you want to select by tapping a drop down button.
- To apply the relevant information, tap the 'Apply' button.
- Proceed with 'Cal 2' following the same steps.




	<ul style="list-style-type: none"> ■ Tap 'QC 1' button, you can select item and lot number that you want to select by tapping a drop down button. ■ To apply the relevant information, tap the 'Apply' button. ■ Proceed with 'QC 2' following the same steps. 	
<p>8.1.12</p>	<ul style="list-style-type: none"> ■ When selection of all calibrators and controls are completed, the relevant information is displayed. 	
<p>8.1.13</p>	<ul style="list-style-type: none"> ■ Tap 'Start' button to start calibration and QC. ■ 'Start' window will pop-up. ■ Check the lot number and the test item. ■ To confirm and proceed with test, tap 'Confirm' button. 	

<p>8.1.14</p>	<ul style="list-style-type: none"> When the test starts, ez-Track¹ check the test item and lot number automatically. The name of the test is displayed and 'In progress' message will be presented below. While the test is in progress, the remaining time is displayed. 	
<p>8.1.15</p>	<ul style="list-style-type: none"> When the test is completed, ez-Track¹ will generate a buzzer signal, and the slot comes out. The test results are saved in the internal memory automatically, and the thermal printer will print out the test results automatically. 	<pre> ----- DATE: Jun-11-2016 / 13:57:52 INSTRUMENT S/N: TPRR020LJ003 USER ID: admin PCT REAGENT LOT: PCMBA13F CALIBRATION DONE CAL F: a(0.82) b(0.13) CAL1: 11860 CAL2: 74400 CAL2: 121800 CAL2: 129300 </pre> <p><Printed Calibration test result></p> <pre> ----- DATE: Jun-11-2016 / 13:59:15 PATIENT ID: PCT 7.58 mg/L SAMPLE: QC1 INSTRUMENT S/N: TPRR020LJ003 USER ID: admin REAGENT LOT: PCMBA13F CALIBRATED_ON: 2016-Jun-11/13:57:52 </pre> <p><Printed QC1 test result></p> <pre> ----- DATE: Jun-11-2016 / 13:58:52 PATIENT ID: PCT 90.81 mg/L SAMPLE: QC2 INSTRUMENT S/N: TPRR020LJ003 USER ID: admin REAGENT LOT: PCMBA13F CALIBRATED_ON: 2016-Jun-11/13:57:52 </pre> <p><Printed QC2 test result></p>


8.1.16

- To check the calibration results, tap the setting icon at the bottom of the Cal/QC screen.
- Tap the 'Cal results' menu. Then ez-Track¹ will display the available test items.
- Tap a test item for which you want to check the results.



8.1.17

- ez-Track¹ displays the lot number of cartridges and calibrators and the next calibration date.
- To check the detailed calibration result, tap the item that you want to check.
- 'Result details' window will pop up and shows the detailed test information as follows.
 - Date / Instrument serial no. / User ID / Name of test item / Reagent Lot / Calibration results
- If you want to print out this information, tap the 'Print' button. Then ez-Track¹ will print out the detailed test information.



Cal results - PCT

Reagent Lot	ST Lot	Next Cal
PCMB13F	PCCALG01_1(Lv.1)	Aug-11-2016

Result details

DATE: Jun-11-2016 / 13:57:52
INSTRUMENT S/N: TPRR020J003
USER ID: admin

PCT

REAGENT LOT: PCMB13F

CALIBRATION DONE

CAL F: a(0.82) b(0.13)
CAL1: 11860
CAL1: 74400
CAL2: 121800
CAL2: 129300

Print Cancel

DATE: Jun-11-2016 / 13:57:52
INSTRUMENT S/N: TPRR020J003
USER ID: admin

PCT


REAGENT LOT: PCMB13F





CALIBRATION DONE




CAL F: a(0.82) b(0.13)
CAL1: 11860
CAL1: 74400
CAL2: 121800
CAL2: 129300

8.1.18

- To check the QC results, tap the setting icon at the bottom of the Cal/QC screen.
- Tap the 'QC results' menu. Then ez-Track¹ will display QC1 and QC2.
- Tap QC1 or QC2. Then ez-Track¹ will display the available test items for QC1 and QC2.
- Tap a test item for which you want to check the results.




<p>8.1.19</p>	<ul style="list-style-type: none"> ■ Checking QC1 ■ ez-Track¹ displays the lot number of cartridges, QC date & time and QC result. ■ If you want to check the QC history, tap  icon. 	 
<p>8.1.20</p>	<ul style="list-style-type: none"> ■ To check the detailed QC1 result, tap the item that you want to check(refer to clause 8.1.18). ■ 'Result details' window will pop up and show the detailed test information as follows. ■ If you want to print out this information, tap the 'Print' button. 	

	<ul style="list-style-type: none"> ■ Then ez-Track¹ will print out the detailed test information. <ul style="list-style-type: none"> - QC Date / Name of test item / Test value / Instrument serial no. / User ID / Reagent Lot / Calibration date 	<div style="border: 1px solid black; padding: 5px;"> <p>-----</p> <p>DATE: Jun-11-2016 / 13:59:15</p> <p>PCT</p> <p style="padding-left: 20px;">7.58 mg/L</p> <p>SAMPLE: QC1</p> <p>INSTRUMENT S/N: TPRR020LJ003 USER ID: admin REAGENT LOT: PCMBA13F CALIBRATED_ON: 2016-Jun-11/13:57:52</p> </div>						
<p>8.1.21</p>	<ul style="list-style-type: none"> ■ Checking QC2 ■ ez-Track¹ displays the lot number of cartridges, QC date & time and QC result. ■ If you want to check the QC history, tap  icon. ■ ez-Track¹ displays the QC2 history as presented on the right. ■ If you want to go back to the previous screen, tap the  icon. 	 <p>The top screenshot shows the 'QC2 - PCT' screen with a table of results:</p> <table border="1"> <thead> <tr> <th>Reagent Lot</th> <th>Date&Time</th> <th>Results</th> </tr> </thead> <tbody> <tr> <td>PCMBA13F</td> <td>Jun-11-2016 / 13:58:52</td> <td>PASS</td> </tr> </tbody> </table> <p>The bottom screenshot shows the 'QC history' screen with a line graph. The y-axis represents concentration in mg/L, ranging from 40 to 160. The x-axis represents time in minutes, ranging from 1 to 12. A single data point is visible at approximately 100 mg/L at the 1-minute mark.</p>	Reagent Lot	Date&Time	Results	PCMBA13F	Jun-11-2016 / 13:58:52	PASS
Reagent Lot	Date&Time	Results						
PCMBA13F	Jun-11-2016 / 13:58:52	PASS						

8.1.22

- To check the detailed QC2 result, tap the item that you want to check(refer to clause 8.1.18).

- 'Result details' window will pop up and show the detailed test information as follows.
- If you want to print out this information, tap the 'Print' button. Then ez-Track¹ will print out the detailed test information:
 - QC Date / Name of test item / Test value / Instrument serial no. / User ID / Reagent Lot / Calibration date



```

DATE: Jun-11-2016 / 13:58:52

PCT
    90.81 mg/L

SAMPLE: QC2

INSTRUMENT S/N: TPRR020LJ003
USER ID: admin
REAGENT LOT: PCMB13F
CALIBRATED_ON: 2016-Jun-11/13:57:52
    
```

8.1.23

- Tap the 'Calibration Alarm' menu. The check box will be toggled whenever you tap this menu. If you want to get an alarm for a past due calibration cycle, the check box should be checked with .


- To change a calibration cycle, tap the 'Calibration cycle' menu.

- ez-Track¹ displays the 'Calibration cycle' window.

- You can change a calibration cycle by tapping the dropdown button in the 'Calibration cycle' window. To confirm the change of a calibration cycle, tap the 'Apply' button.

Note: You can change a calibration cycle by day, week and month.


- If a calibration cycle is past due, ez-Track¹ displays 'Alarm details' window. Also a red colored 'A' marking is displayed on the setting icon under the main screen. In this case, you need to perform a new calibration.





8.1.24

- 'Enable Report' is for sending Cal/QC data to the server.
- Tap the 'Enable Report' menu to activate this function.
- ez-Track¹ displays 'Enable Report' window for confirmation of sending Cal/QC data to Boditech's server only for maintenance purposes.
- If you agree with the use of this function, tap the 'Agree' button.
- If you do not agree, tap the 'Cancel' button.
- The check box will be checked as .


Note: When the 'Enable Report' menu is activated and ez-Track¹ is connected to LAN, ez-Track¹ will send Cal/QC data and the device usage history to Boditech's server automatically. This information is used only for maintenance purposes by the manufacturer.



8.2 Retrieval/Printing of stored test results

<p>8.2.1</p>	<ul style="list-style-type: none"> ez-Track¹ can save up to 5,000 test results in its internal memory. Note: If the number of stored test results exceeds 5,000 results, the oldest result will be deleted when a new test result is saved. To print out stored test result(s), tap 'History'. Note: Stored test results cannot be printed if a test is in progress. 	
<p>8.2.2</p>	<ul style="list-style-type: none"> Upon tapping "History", datewise log of stored test results is displayed in a reverse chronological order. Select the test result(s) you want to print out. It is possible to select up to 5,000 test results. Tap 'Print' to print out the selected test result(s). A new window will pop up. Tap "OK" to confirm printing action. If you do not intend to print out the selected test results, tap 'Cancel'. 	

- After tapping 'OK' button, ez-Track¹ will display the printing status and print out selected result(s).



8.2.3

- If you want to search specific test result(s), tap "Search" and input the search term (e.g. date, patient ID, user name, test name etc.)
- Tap the search icon
- Searched results will be displayed.




8.2.4

- If you just want to note the details of stored test result(s) without printing, tap the field displaying the test result of your interest.
- A new window will pop up displaying details such as patient ID, test name, test result, date & time of the test, lot number and expiry date of test reagents, etc.
- If you want to close this window, tap the 'OK' button.




8.2.5

- Printout generated by the built-in thermal printer contains the test result alongwith relevant information as shown herewith.
- The printout will also contain test-specific reference values/interval for the specific patient population if such information has been previously inputted by the laboratory.



8.2.6


- You can find Date, PID, Item and Result in this History screen.
- If you want to align and check the results by items, please tap the arrow button next to Item.



8.3 Holder In/Out

8.3.1

- The cartridge holder will go inside or come out by tapping the Holder In/Out button at the bottom of the screen.
- The In/Out button cannot be used while a test is in progress.



9. Configuring

Following settings and/or parameters of ez-Track¹ can be adjusted.



Clause 9.1 and 9.2 are only accessible to the 'Administrator' only.

9.1 Configuring System Operation



- 9.1.1 Date, Time, Date & Time format
- 9.1.2 Language
- 9.1.3 Printer Option
- 9.1.4 Test menu
- 9.1.5 Device Information
- 9.1.6 Test counter
- 9.1.7 Setting Tip / Sample type
- 9.1.8 LIS Configuration
- 9.1.9 Updating the S/W



9.2 User management



9.3 ID chip management



9.4 Temperature



9.1 Configuring System Operation



<p>9.1.1.1</p>	<p>9.1.1 Date, Time Date & Time format</p> <ul style="list-style-type: none">■ Setting the 'Date':<ul style="list-style-type: none">① Tap the 'Configuration menu'.② Tap the 'Set date' menu.	 <p>The screenshot shows the 'System' configuration screen. At the top, there is a blue header with a back arrow and the word 'System'. Below it, the 'Configuration menu' is highlighted with a red dashed box. Underneath, there are 'Miscellaneous functions' and a large grey rectangular area. Below that, the 'Configuration' screen is shown, with the 'Set Date' option highlighted by a red dashed box. The 'Configuration' screen lists several settings: Set Date, Set Time, Set Date & Time Format, Ethernet Configuration, LIS Configuration, Language - English (United States), Printer Option, and Set Default Type.</p>
<p>9.1.1.2</p>	<ul style="list-style-type: none">■ Tap the + or - button(s) to set the desired month, date and year.■ To set the changed date, tap 'Set'.	 <p>The screenshot shows the 'Set date' dialog box. It has a dark background with a title bar that says 'Set date'. Below the title bar, there are three columns of controls. The first column has a '+' button above 'Dec' and a '-' button below. The second column has a '+' button above '06' and a '-' button below. The third column has a '+' button above '2016' and a '-' button below. At the bottom, there are 'Set' and 'Cancel' buttons. The '+' buttons for month, date, and year are highlighted with red dashed boxes.</p>



<p>9.1.1.3</p>	<ul style="list-style-type: none">■ Setting the 'Time':<ol style="list-style-type: none">① Tap the 'Configuration menu'. ② Tap the 'Set time' menu.	 <p>The screenshot shows the 'System' menu with a blue header and a back arrow. Below the header, 'Configuration menu' is highlighted with a red dashed box. Underneath, there is a section for 'Miscellaneous functions' which is currently greyed out.</p>
<p>9.1.1.4</p>	<ul style="list-style-type: none">■ Tap the + or - button(s) to set the desired values for hour and minute.■ You can set the time as 'AM' or 'PM' by tapping the adjacent tab.■ To set the changed time, tap 'Set'.	 <p>The screenshot shows the 'Configuration' menu with 'Set Time' selected. A 'Set time' dialog box is open, showing a digital clock interface. The hour is set to 2 and the minute to 25. The 'PM' tab is selected. The '+' and '-' buttons for both hours and minutes are highlighted with red dashed boxes. At the bottom of the dialog, the 'Set' button is also highlighted with a red dashed box.</p>



<p>9.1.1.5</p>	<ul style="list-style-type: none">■ Setting the Date Format① Tap the 'Configuration menu'. ② Tap the 'Set Date & Time format' menu.	
<p>9.1.1.6</p>	<ul style="list-style-type: none">■ Check/tap the blank box corresponding to the date format you want to select.■ If you want to cancel changing the date format, tap 'Cancel'.	




<p>9.1.2.1</p>	<p>9.1.2 Language</p> <ul style="list-style-type: none"> ■ Setting the Language: <p>① Tap the 'Configuration menu'.</p> <p>② Tap the 'Language' menu.</p>	
<p>9.1.2.2</p>	<ul style="list-style-type: none"> ■ Select the language you want to use. ■ The system reboot process is required for the selected language to be applied. 	



<p>9.1.3.1</p>	<p>9.1.3 Printer Option</p> <ul style="list-style-type: none"> ■ Turning the Printer On / Off: <p>① Tap the 'Configuration menu'.</p> <p>② Tap the 'Printer Option' menu.</p>	
<p>9.1.3.2</p>	<ul style="list-style-type: none"> ■ Check/tap the blank box of the 'Printer on' menu to turn the printer 'on'. ■ Uncheck/tap the checked box of the 'Printer on' menu to turn the printer 'off'. 	



<p>9.1.3.3</p>	<ul style="list-style-type: none"> ■ If you want to print out a test result with a 'Reference Value', check/tap the blank box of 'Print including the reference value'. 	
<p>9.1.3.4</p>	<ul style="list-style-type: none"> ■ If you want to edit 'Reference value', tap the 'Edit reference value' menu. ■ When you tap a blank area next to the item name, ez-Track¹ displays an edit window for a reference value and unit. 	


<p>9.1.4 Test Menu</p> <ul style="list-style-type: none"> ■ Setting 'Test Menu': ■ You can test functions for LCD display, LCD touch, LAN connection (only for manufacturer), SD card connection (only for manufacturer) and printer connection. <p>① Tap the 'Miscellaneous functions' menu.</p> <p>9.1.4.1</p> <p>② Tap the 'Test menu'.</p>		
<p>9.1.4.2</p> <ul style="list-style-type: none"> ■ To test LCD display, tap the 'LCD display' menu. ■ ez-Track¹ displays various colors on the screen during LCD display testing. <p>■ If there is no problem with the LCD display, tap the 'Exit' button in the center of the screen. If you find some problems with the LCD display, please contact Theradiag's Client Support or a regional distributor.</p>		

<p>9.1.4.3</p>	<ul style="list-style-type: none"> ■ To test 'LCD touch', tap the 'LCD touch' menu. ■ ez-Track¹ will display 9 blocks(8 'Hide' blocks and 1 'Exit' block). ■ To test 'LCD touch', tap the 8 'Hide' blocks. ■ To complete this test, tap the 'Exit' block. 	
<p>9.1.4.4</p>	<ul style="list-style-type: none"> ■ To test the network connectivity, tap the 'Network' menu. ■ If the LAN is connected, the Network type is displayed as 'Ethernet' and Network status is displayed as 'Connected'. 	

	<ul style="list-style-type: none"> ■ If the LAN is not connected, Network type is displayed as 'None' and Network status is displayed as 'Disconnected'. 	
<p>9.1.4.5</p>	<ul style="list-style-type: none"> ■ To check a SD card, tap the 'SD card' menu. ■ ez-Track¹ will check the SD card inserted in the 'SD card slot'. ■ If SD card has been inserted, screen will display the alert, 'SD card mounted'. ■ If SD card has not been inserted or if the SD card has been inserted yet could not be accessed by the system, the screen will display 'No SD card'. 	
<p>9.1.4.6</p>	<ul style="list-style-type: none"> ■ To test printer connection, tap the 'Printer' menu. ■ ez-Track¹ prints out standard texts according to language setting as shown in the pictures herewith. 	 <pre> PRINT-TEST-START-START-START-*** ABCDEFGHIJKLNPQRSTUWXYZ\ !@.# ABCDEFGHIJKLNPQRSTUWXYZ\ !@.# BCDEFGHIJKLNPQRSTUWXYZ\ !@.# CDEFGHIJKLNPQRSTUWXYZ\ !@.# DEFGHIJKLNPQRSTUWXYZ\ !@.# EFGHIJKLNPQRSTUWXYZ\ !@.# FGHIJKLNPQRSTUWXYZ\ !@.# GHIJKLNPQRSTUWXYZ\ !@.# HIJKLNPQRSTUWXYZ\ !@.# IJKLNPQRSTUWXYZ\ !@.# JKLNPQRSTUWXYZ\ !@.# KLNPNQRSTUWXYZ\ !@.# LNPQRSTUWXYZ\ !@.# MNPQRSTUWXYZ\ !@.# NPQRSTUWXYZ\ !@.# OPQRSTUWXYZ\ !@.# PRINT-TEST-END-END-END-END-*** </pre> <p>English format</p>

	<p>③ Following information related to the system will be displayed:</p> <ul style="list-style-type: none"> ▪ Build Date ▪ Build ID ▪ Android Version ▪ Kernel Version ▪ Application Version ▪ Firmware Version 	
<p>9.1.6</p>	<p>■ Setting the Test counter:</p> <ul style="list-style-type: none"> ■ You can use this function to find out the cumulative number of tests performed in ez-Track¹. ■ You can identify the following: <ul style="list-style-type: none"> - C1 : Total number of tests performed since the initial use - C2 : Total number of tests performed after the latest repair or servicing of the system <ol style="list-style-type: none"> 1. Tap the 'Test counter' menu. 2. To check C1 and C2, tap the 'Reset' menu. <p>■ If you want to reset C2, tap the 'Reset' button.</p>	

<p>9.1.7</p>	<p>■ Setting Default Type</p> <p>This is for setting default Tip/Sample type.</p> <p>In following cases, ez-Track¹ sets Tip/Sample type as the default setting automatically.</p> <ol style="list-style-type: none"> 1. When ez-Track¹ turns on. 2. After finishing a test. 3. When user stops a current test by tapping the 'Stop' button. 	
<p>9.1.8</p>	<p>■ LIS Configuration</p> <p>■ For sending test results from ez-Track¹ to the LIS server via Local Area Network or RS232C.</p> <ol style="list-style-type: none"> 1. Tap the 'LIS configuration' menu. 2. To send test results by LIS after finishing tests, select 'On' in the 'Send' item. In the opposite case, select "Off". 3. Please select the data format. 4. Please select the type of interface you wish to use. <p><i>Note: Once 'Network' is selected as a preferred interface, only HL7 protocol is available and the protocol item will be changed to 'HL7' automatically. Server field is only available when 'Network' is selected. You can input a server IP or PORT in the server field.</i></p> <p>■ To apply setting, please tap the 'Save' button.</p>	


<p>9.1.9.1</p>	<ul style="list-style-type: none">■ Updating Software: <i>Note: This function should be handled only by the 'administrator', a trained engineer or manufacturer's authorized representative. Incorrect or unwarranted use of this function can cause serious problems in system operations.</i>■ If you insert the 'Upgrade SD card' and tap the 'SD/USB' menu, then the program will automatically determine if an upgrade is needed. When an upgrade is needed, it will start upgrading.■ If you want to update S/W via network, you can tap 'Network' menu. Please refer to clause 9.1.9.4	 <p>System</p> <p>Configuration menu</p> <p>Miscellaneous functions</p> <p>Miscellaneous</p> <p>SW update</p> <p>Test menu</p> <p>Device information</p> <p>Bay counter</p> <p>SW update</p> <p>SD/USB</p> <p>Network</p>
----------------	--	--



9.1.9.2







■ **Updating Firmware via a SD card:**

1. Tap the 'Miscellaneous functions' menu.
 2. Tap the 'SW update' menu.
 3. Tap the 'SD/USB' menu.
 4. 'Update via SD-Card' window will pop up.
 5. Tap the 'Update' menu after confirming the information in the pop-up box.
- ① If you insert a SD card with an updated version of firmware, it will start to update automatically.
- ② If you tap the 'Confirm' menu, it will automatically reboot.
- ③ Update will be completed after rebooting.

Note: You can check whether the firmware update is successful by accessing the 'Device information' menu as described in clause 9.1.5.



		
<p>9.1.9.3</p>	<p>■ Updating the Android Operating System via SD card:</p> <p>To update the 'Android OS', insert the 'Upgrade SD Card' (provided by the manufacturer or its authorized service representative) into the 'SD card port'.</p> <ol style="list-style-type: none"> 1. Tap the 'Miscellaneous functions' menu. 2. Tap the 'SW update' menu. 3. Tap the 'SD/USB' menu. 4. 'Update via SD-Card' winddjow will pop up. 5. Tap the 'Update' menu after confirming the information in the pop-up box. <ol style="list-style-type: none"> ① If you insert a SD card with an updated version of Android OS, it will start to update automaticallly. ② If you tap the 'Confirm' menu, it will automatically reboot. ③ Update will be completed after rebooting. 	

	<p>6. Then "SW update" screen will be presented. To finish Android OS update, it takes about 4 minutes.</p> <p>7. When Android OS update is completed, ez-Track¹ will reboot automatically.</p> <p><i>Note: You can check if the Android OS update is successful by accessing the 'Device information' menu as described in clause 9.1.5.</i></p>	    
<p>9.1.9.4</p>	<p>■ Updating the S/W via Network:</p> <p>To update the S/W via network, ez-Track¹ should be connected to network (e.g. LAN cable connection is needed.) and insert a 'SD Card' (provided by the manufacturer or its authorized service representative) into the 'SD card port'.</p>	








1. Tap the 'Network' menu.
2. If the ez-Track¹ is not registered, 'Device Registration' window will be displayed.
3. Tap the 'Register now' button to register the device.

Note: If you tap the 'Later' button, It will not update S/W via network.



4. Input registration information and tap the 'OK' button. Then tap the 'Network' menu again(refer to Step 1).
5. ez-Track¹ displays 'Update via Network' window and you can check the available updates.
6. Tap the 'Download' button to proceed with the update.
7. ez-Track¹ displays 'Download update' window and the downloading progress is displayed.
8. When the download is completed, ez-Track¹ performs the update automatically.



9.2 User management

<p>9.2.1</p> <ul style="list-style-type: none"> ■ 'User Management' is accessible to the 'Administrator' only. ■ This function should be used for adding a new user or deleting an existing/registered user. ■ The 'Administrator' can register/add up to 100 users at a given time. <ol style="list-style-type: none"> 1. Tap the 'User management' menu. <i>(If any registered user (xxxx) other than the 'administrator' tries to access the 'User management' menu, ez-Track¹ will display the warning message 'Access denied for xxxx'.)</i> 2. Tap the '+' icon to add a user. 3. Input suitable 'user name' and 'password'. 4. Confirm/re-enter the password. <p> <i>User name and password should comprise of at least 4 characters each.</i></p> <ol style="list-style-type: none"> 5. Tap 'OK' to add the user with above user name and password. 6. Tap 'Cancel' if you do not intend to proceed. 	  
<p>9.2.2</p> <ul style="list-style-type: none"> ■ To delete an existing/registered user, check/select the blank box corresponding to the 'user name' you want to delete. ■ Tap 'Delete'. ■ Tap 'OK' to confirm deleting the selected user. <p> <i>'Admin' being the default user name, cannot be deleted or modified.</i></p>	 


9.3 ID Chip Management

<p>9.3.1</p>	<ul style="list-style-type: none"> ■ When any 'ez-Tracker Test ID Chip' is inserted into the ID chip slot, relevant information in the ID chip will be saved in the internal memory of ez-Track¹. Note: ez-Track¹ can save ID chip information from up to 500 ID chips. ■ By accessing the 'ID chip management' menu, you can search and check relevant ID chip information such as name of ez-Tracker test, Lot number, Expiry date and sample type(s) ■ You can also delete an ID chip information that is now irrelevant or outdated on account of inventory with the particular lot number being exhausted. ■ Tap the 'ID chip management' menu. 	 <p>The screenshot shows a mobile application interface with a top navigation bar containing 'Analyses', 'History', and 'Cal/QC'. Below this is a 'System' menu with options: 'User management', 'ID chip management' (highlighted with a red dashed box), 'Temperature', and 'Maintenance'. At the bottom, there is a status bar with 'admin', a settings icon, a battery icon, and the date/time 'Jun-12-2017 / 16:00'.</p>								
<p>9.3.2</p>	<ul style="list-style-type: none"> ■ To search ID chip information, tap the 'Search' icon. ■ Input appropriate term/keyword (e.g. test name, lot number, sample type etc.) for the information you want to search and tap 'search'. 	 <p>The screenshot shows the 'ID chip management' screen with a search bar at the top containing the text 'Search'. Below the search bar is a table with 3 results found. The table has two columns: 'Item' and 'Lot No'. The rows are: CRP (CRLCA05M), PCT (PCIXA07M), and HBsAg (ABYKA06). A red dashed box highlights the search bar.</p> <table border="1" data-bbox="732 949 1010 1045"> <thead> <tr> <th>Item</th> <th>Lot No</th> </tr> </thead> <tbody> <tr> <td>CRP</td> <td>CRLCA05M</td> </tr> <tr> <td>PCT</td> <td>PCIXA07M</td> </tr> <tr> <td>HBsAg</td> <td>ABYKA06</td> </tr> </tbody> </table>	Item	Lot No	CRP	CRLCA05M	PCT	PCIXA07M	HBsAg	ABYKA06
Item	Lot No									
CRP	CRLCA05M									
PCT	PCIXA07M									
HBsAg	ABYKA06									


9.3.3


- List of the search results will be displayed.
- To get the ID chip information, tap the name of the test listed under the column 'Item'.

- Relevant information (test name/item, lot number, expiry date and sample type) corresponding to the ID Chip will be displayed.



9.3.4

- To delete the ID chip information, select/check the blank box corresponding to the ID chip.
- Tap the  icon.
- Tap 'OK' to confirm deleting the ID chip information.



9.4 Temperature


- You can access this function to check the operating temperature.

Note: The normal/operational temperature range is 33~37 °C.

1. Tap the 'Temperature' menu.
2. Operating temperature status of ez-Track¹ will be displayed.

9.4.1

Note: If operating temperature status is being shown as 'Incorrect temperature' or 'Error status', you cannot perform ez-Tracker tests on the system until the ambient temperature is managed to a value in the operational range.




The screenshot shows the 'Maintenance' menu with 'Temperature' selected. The top navigation bar includes 'Analyses', 'History', and 'Cal/QC'. The user is logged in as 'admin' on July 12, 2017, at 18:00. The device ID is AA1.WW.02.05.03.

Units	Temperature	Status
°C	34.8	●
°F	94.6	●

Legend:
● - Ok, no errors
! - Incorrect temperature
● - Error status

<Normal status>



The second screenshot shows the 'Temperature' screen with an error status. The temperature is 30.5 °C (86.9 °F), and the status is indicated by an orange exclamation mark icon.

Units	Temperature	Status
°C	30.5	!
°F	86.9	!

Legend:
● - Ok, no errors
! - Incorrect temperature
● - Error status

<Incorrect status>

10. Maintenance, Servicing and Disposal


10.1 Maintenance



No maintenance other than the following is required for ez-Track¹:

- Loading paper in the built-in thermal printer
- Emptying the 'Waste Tip Box' and resetting the 'Waste tip counter'
- Occasional cleaning of cartridge holders and exterior
- Raw data backup for maintenance or service purposes



Loading Paper in the Printer:

- Open the printer cover as shown in the figure below, insert the paper in the proper direction, and take out the paper at least 1 cm.
- Remove the remaining paper, if any.



	When the printer cover is closed, built-in printer ejects the printer paper a little for initialization.
	Be careful exerting excessive force when opening or closing the printer cover for preventing malfunction of built-in printer.


Emptying the 'Waste Tip Box' and Resetting the waste tip counter:

- To reset the counter indicating the collection status of wasted pipette tips, pull out the waste tip box.
- The 'Reset waste box' will pop-up.
- Empty the 'Waste Tip Box' according to the instructions in the pop-up window.
- To reset the waste tip counter, tap  button.
- If you want to terminate this procedure without a reset, tap  button.
- If reset is completed, the percentage under the trash can icon will be presented as "0 %".




Warning message when user does not empty the 'Waste box':


- ez-Track¹ displays a warning pop-up window for recommending that user needs to empty the waste tip box before starting a test. In this case, you need to empty the waste tip box and reset the waste tip counter by tapping the icon.
- If you proceed with a test without emptying the waste tip box and resetting the waste tip counter, ez-Track¹ will lock all functions for test when the waste tip counter reaches to the 120%.




- **Raw data backup:**
- Raw data for the tests are needed to identify the cause of problems and solve them when these problems occur.
- Tap the **'Extension menu'** in the main menu.




- Tap the **'Maintenance'** menu.



- Tap the 'Raw data' menu.




- Insert a SD card into the SD card port.
- Select raw data that you want to backup by tapping the check box.
- Tap the 'Backup' button.




- Tap the 'OK' button to proceed with backup.

Note: *The backup data is saved in boditech/raw data folder. This data is encoded and it is only accessible by the manufacturer.*



10.2 Cleaning

- Occasionally clean the cartridge holders and exterior of ez-Track¹ by wiping with a soft cloth slightly dampened with 70% alcohol or 0.5% bleaching solution.
- Do not use soap solution or any industrial solvent for cleaning purpose.


	Unplug ez-Track ¹ from external power source before cleaning.
---	--

10.3 Servicing

- ez-Track¹ does not have any operator-serviceable component.
- If any technical servicing or maintenance is required, ez-Track¹ should be sent back to Theradiag or its authorized service agent.



10.4 Disposal

- If ez-Track¹ is to be disposed of for any reason, the user is advised to follow local regulations as applicable to the disposal of electrical equipment.

	At the end of its useful life, ez-Track ¹ could be sent back to Theradiag for recycling or proper disposal. Alternatively, it should be disposed of in accordance with local regulations.
---	---

10.5 Transportation and Storage

- The original shipping container/box should be used to ship or transport ez-Track¹. The same container/box is also recommended for an extended storage period.

	When transporting or storing ez-Track ¹ , keep it dry in an upright position and protect it from mechanical shocks.
	Before transporting ez-Track ¹ , remove tips from the tip box.

11. Trouble shooting

ez-Track¹ may encounter some operational problems occasionally. The symptoms/error codes displayed, probable causes and recommended remedial/corrective measures for these problems have been listed in the table below.

Symptoms	Probable Causes	Recommended Remedial/ Corrective Measures
Nothing happens; ez-Track ¹ does not show any sign of functioning	Power failure	Check whether external power outlet is live.
	Poor connection between the power adapter and the power cable.	Remove the power cable and reconnect firmly.
	Faulty power adapter	Call Theradiag's Technical Services.
	Faulty power cable	Call Theradiag's Technical Services.
	Main Power Switch is 'Off'.	Turn the main power switch 'On'.
	ez-Track ¹ is in 'sleep' mode	Gently press the power button located on the front of ez-Track ¹ . In the problem still persists, call Theradiag's Technical Services.
'Cartridge Holder' does not protrude out as and when required.	Mechanical failure	Call Theradiag's Technical Services.
	Software malfunction	Turn the power 'Off' and turn it 'On' again.
'LCD Display' does not function properly.	Electrostatic shock	Connect to a grounded external power outlet.
		Remove any source of electrostatic charge. Reset the system.
Clicking noise generated while the cartridge holder(s) protrudes out	Electrical circuit failure	Call Theradiag's Technical Services.
	Normal mechanical setting.	No corrective measures are not required.

Barcode error	ez-Track ¹ can't read the barcode of the ez-Tracker test cartridge on account of the barcode being misprinted or obscured due to presence of deformation or foreign matter.	Remove any foreign matter obscuring the barcode and continue the test by reinserting the same cartridge. If the error code persists, repeat the test using a new 'ez-Tracker test cartridge'.
	Malfunction of the barcode- reading/scanning module.	Call Theradiag's Technical Services.
Insert printer paper	No paper in the built-in printer.	Load a new paper roll.
Error 25-29	Technical problem with the 'ez-Tracker test Cartridge' or malfunction.	Repeat the test using a new 'ez-Tracker test Cartridge'.
Insert tip	Test run has been initiated without placing a 'pipette tip' into the 'tip slot' of the ez-Tracker test cartridge inserted in the cartridge holder.	Place a 'pipette tip' into the 'tip slot' of the sample-loaded 'ez-Tracker test cartridge' (inserted in the cartridge holder) before initiating the test run (i.e. before tapping 'Start' on the screen display).
	ez-Track ¹ can't detect the 'pipette tip' placed into the 'tip slot' of the inserted test cartridge.	Call Theradiag's Technical Services.
Please wait until the reader reaches recommended temperature.	System temperature is outside the permissible or operational temperature range.	After pressing the 'stop' button, check the temperature status of the system (Please refer to 'Clause 9.4: Temperature'.) If system temperature is outside the operational range, wait for about 15 minutes. Proceed with the test(s) after

		temperature status of the system changes to 'normal'.
ERR_60	Malfunction of the built-in temperature sensor.	Call Theradiag's Technical Services.
Insert ID chip	System check cartridge is not inserted.	Insert a system check cartridge and do a system check.
	The test is performed without inserting ID chip.	Insert an ID chip and do the test again.
	Lot number of ID chip and test cartridge do not match.	Perform the test again with matching ID chip and test cartridge.
Error_S1	System check is not successful.	Try to perform another round of system check.
	Problems in the optical parts.	If this error is repeated, call Theradiag's Technical Services.
The password does not match.		Input correct password (for administrator).
	Inputed password does not match with the administrator's password.	If administrator has forgotten his/her password, call Theradiag's Technical Services for resetting administrator's password.
	Inputed password does not match with the registered user's password.	Input correct password or if the user has forgotten his/her password, reset the password with the help of 'Administrator'.
Built-in printer does not print the test results	The printer has not been set 'on'	Set the printer 'on' by accessing the 'Printer on/off' menu. (Refer to Clause 9.1.9~9.1.10).
	Malfunction of the printer.	Call Theradiag's Technical Services.

12. Warranty

ez-Track¹ has passed strict quality assurance and testing procedures.

Theradiag's expressed and implied warranties are conditional upon full observance of Theradiag's published directions with respect to the use of Theradiag's products.

To avail the warranty service, the consumer must return the defective ez-Track¹ unit or any component thereof; along with proof of purchase to Theradiag or its nearest Authorized Service Agent.

Warranty Information

During the warranty period (24 months from the date of purchase), defective or spontaneously malfunctioning ez-Track¹ or its component(s) will be repaired free of cost or may be compensated for in accordance with consumer protection rules and regulations.

The consumer may be charged for repair of ez-Track¹ even during the warranty period in following cases:

- Improper use or misuse
- Consumer's intentional mishandling, abuse or neglect of the product(s)
- Unauthorized repair / replacement of or alteration in any component
- Damages or defects due to fire, pollution, earthquake or other natural disasters.

For after-sales inquiry, please contact:

Theradiag

14 Rue Ambroise Croizat 77183 Croissy Beaubourg France

Tel: +(33)-1-64-62-10-12 / Fax : +(33)-1-64-62-09-66

e-mail: sav@theradiag.com

Website: www.theradiag.com

Warranty Card

Thank you for purchasing ez-Track¹.

Please fill the required information and send a copy of this warranty card within 30 days of date of purchase to the above address.

Product Name	ez-Track ¹
Catalog No.	ET 001
Serial No.	
Date of Manufacture	
Warranty Period	24 months from the date of purchase
Date of Purchase	Year _____ Month _____ Day _____
Manufacturer	Boditech Med Inc.
Distributor	Theradiag
Purchase Location	

ez-Track¹



Boditech Med Inc.

43, Geodudanji 1-gil, Dongnae-myeon,
Chuncheon-si, Gang-won-do 24398, Korea
Tel: +82-33-243-1400 Fax: +82-33-243-9373
강원도 춘천시 동내면 거두단지 1길 43



Obelis s.a.

Bd. Général Wahis 53, 1030 Brussels, Belgium
Tel: +32-2-732-59-54 Fax: +32-2-732-60-03
mail@obelis.net



Theradiag

14 Rue Ambroise Croizat 77183 Croissy Beaubourg France
Tel: +33-1-64-62-10-12 Fax: +32-1-64-62-09-66
www.theradiag.com