

**Updated!**  
CardioChek Plus  
Analyzer  
REF 2700

# Lipid Panel and Glucose results in as little as 90 seconds

## Now with:

- Updated interface
- Quick backout options
- Simplified navigation
- Test counter



**Lipid Panel**  
Test Strip  
REF 1710

**eGLU**  
Test Strip  
REF 2713

## Lab-quality results at the point of care

### Accurate

Total and HDL Cholesterol meet NCEP guidelines for accuracy and precision

### Portable

Battery-operated  
Refrigeration not required for test strips or controls

### User-friendly

Fingerstick blood sample  
CLIA-waived  
Self-guided online training available

### Cost-effective

Minimal start-up costs  
No standard analyzer contracts  
Tests may be reimbursable

### Measures

- Total cholesterol
- HDL cholesterol
- Triglycerides
- Glucose

### Calculates

- LDL cholesterol
- TC/HDL ratio
- LDL/HDL ratio
- Non-HDL cholesterol

# The PTS Panels Lipid Panel+eGLU Smart Bundle 2-Pack

For simultaneous lipid panel and glucose results



## REF 2729

Includes 15 lipid panel test strips; 15 eGLU test strips

	REF#	Analytes
<b>Three-Analyte</b> Use PTS Collect™ 40 µL capillary tubes (REF 2866) with three-analyte test strips.	1710	● Lipid Panel
	2412	● CHOL+HDL+GLU
<b>Dual-Analyte</b> Use PTS Collect 30 µL capillary tubes (REF 2865) with dual-analyte test strips.	1765	● CHOL+GLU
	1821	● CHOL+HDL
<b>Single-Analyte</b> Use PTS Collect 15 µL capillary tubes (REF 2863) with single-analyte	1711	● Cholesterol
	1716	● Triglycerides
	1713	● Glucose
<b>No capillary tube required</b> with eGLU test strips.	2713	● eGLU

## Controls and Connectivity

### PTS Panels® Multi-Chemistry Controls

REF 721



- Ready-to-use control materials
- Two-level control kit
- Room temperature storage\*

\* Unopened controls can be stored at room temperature until expiration date. Opened controls can be stored at room temperature for 10 months, or until the expiration date, whichever comes first.

### PTS Connect® Printer

REF 2750



- Outputs results on self-adhesive labels
- Can connect to CardioChek Plus analyzer through a provided cable.



For more information or to place an order:

1 (317) 870-5610  
 customerservice@ptsdiagnostics.com  
 ptsdiagnostics.com