

BSX READER® INSTRUMENTS



The BSX READER® provides standardised and automated interpretation of lateral flow tests from the BSX range

🕒 ***Easy to use and intuitive***

The touch screen guides you through each step of the result interpretation

🕒 ***Full traceability***

Operator, patient and cassette batch traceability allowing you to meet your quality objectives with peace of mind

🕒 ***Standardized reading of results***

Automated reading and interpretation eliminating user error in the interpretation of results

BSX READER®

The BSX READER® is a highly **sensitive and robust** measuring system for the lateral flow tests provided by BIOSYNEX. This **flexible and reliable** system allows users to easily run the lateral flow tests of the **BSX Range**.

The instrument will guide you through each step of the result interpretation by means of a series of messages and animations displayed on the touch screen.

BSX READER®

Easy to use and intuitive operation

Thanks to the touch screen

Adaptable to multiple cassette formats

Thanks to the integrated adaptive drawer

Automatic or manual incubation

For each test in the BSX Range

Ergonomic design

Compact and portable reader



Auto-calibration

During the start-up

Wide range of compatible tests

Dedicated BSX range

Optimization of technician time

No maintenance

Connectable to LIS

HL7 Protocol

BSX Range

Women's health:

AMNIOQUICK® Card BSX
AMNIOQUICK® DUO+ BSX
PREMAQUICK® BSX

Bacteriology:

TETANOTOP® BSX
BIOSYNEX S. pneumoniae BSX
BIOSYNEX L. pneumophila BSX
BIOSYNEX STREP A BSX

Virology:

BIOSYNEX® COVID-19 Ag BSX
INFLUENZATOP® BSX
RSVTOP® BSX

Biochemistry:

DRUGCHECK®+ 11 BSX
DRUGCHECK®+ 10 BSX
DRUGCHECK®+ 5 BSX
DRUGCHECK®+ THC BSX

The BSX READER® is intended to be used only in combination with the lateral flow tests of the BSX range and only for the applications described by BIOSYNEX. For more information, please contact your local BIOSYNEX sales representative.

BSX READER® Kit content

REF. 5060029

- BSX READER® instrument
- Power supply cable
- Barcode reader
- User manual

ACCESSORIES

- Dymo LabelWriter 450 printer (ref. 502250)
- Printing labels

Manufacturer: BIOSYNEX S.A. (France)

In vitro diagnostic medical devices.

Read the instructions in the data sheets carefully.

FLY_5060029_EN_V05202203R01
Rev : 03/22

BIOSYNEX

BIOSYNEX S.A.
22 Boulevard Sébastien Brant
67400 ILLKIRCH-GRAFFENSTADEN - France
Phone : +33 3 88 78 80 88 - Fax : 03 88 78 76 78
E-mail : export@biosynex.com
www.biosynex.com