




BIOSYNEX

RAPID
DIAGNOSTIC
TESTS

WWW.BIOSYNEX.COM

A comprehensive range of products for the diagnosis of seasonal respiratory infections

BIOSYNEX® COMBO COVID-FLU BSS

QUALITATIVE DETECTION OF SARS-COV-2, INFLUENZA A & INFLUENZA B VIRUS ANTIGENS

- ✓ **FAST:**
Result in 10 minutes
- ✓ **USER-FRIENDLY:**
Test kit with 25 tests containing all necessary materials to perform the test
- ✓ **SAMPLE COLLECTION:**
Nasopharyngeal swab
- ✓ Detection of **Influenza virus Type A** (subtypes **H1N1**, **H3N3**) & **Influenza virus Type B**

BIOSYNEX® COMBO COVID-FLU BSS. Manufacturer: BIOSYNEX SWISS SA (Switzerland).
Packaging: 25 tests + 25 swabs + 25 pre-filled extraction tubes + 25 dropper caps.
Positive controls swabs included for SARS-CoV-2, Influenza A and Influenza B.
Storage: 2 - 30 °C.
For further information, please refer to the instructions for use of the product.



REF : SW40008

BIOSYNEX® COMBO COVID-FLU-RSV BSS

MULTIPLEX TEST FOR QUALITATIVE DETECTION AND DIFFERENTIATION OF SARS-COV-2, INFLUENZA VIRUS (A&B) & RESPIRATORY

- ✓ **FAST:**
Result in 10 minutes
- ✓ **USER-FRIENDLY:**
Test kit with 25 tests containing all necessary materials to perform the test
- ✓ **SAMPLE COLLECTION:**
Nasopharyngeal swab

BIOSYNEX® COMBO COVID-FLU-RSV BSS. Manufacturer: BIOSYNEX SWISS SA (Switzerland).
Packaging: 25 tests + 25 swabs + 25 pre-filled extraction tubes + 25 dropper caps.
Positive controls swabs included for SARS-CoV-2, Influenza A, Influenza B and RSV.
Storage: 2 - 30 °C.
For further information, please refer to the instructions for use of the product.



REF : SW400011



ADVANTAGES

- Determine the cause of respiratory symptoms
- Take appropriate health protection measures
- Manage child protection in case of positive RSV diagnosis

BIOSYNEX® STREP A

+QUALITATIVE DETECTION OF GROUP A BETA-HEAMOLYTIC STREPTOCOCCI

- ✓ **FAST:**
Result in 5 minutes
- ✓ **USER-FRIENDLY:**
Simple procedure in 4 steps
- ✓ **SAMPLE COLLECTION:**
Oropharyngeal swab

BIOSYNEX® Strep A. Manufacturer: BIOSYNEX SWISS SA (Switzerland).
Packaging: 20 tests + 20 swabs + 20 tongue depressors + 2 extraction swabs + 20 extraction tubes.
Positive & negative controls included
Storage: 2 - 30 °C.
For further information, please refer to the instructions for use of the product.



REF : 1030005



ADVANTAGES

- STREP A: Main cause of upper respiratory tract infections
- Differentiation between bacterial and viral angina
- Optimization of antibiotic prescription

Other products related to respiratory infections

BIOSYNEX® COVID-19 Ag BSS




- User-friendly & easy-to-use
- Result in **15 min.**

Evaluated
by the Paul
Ehrlich
Institute

Ref. SW40006



BIOSYNEX® INFLUENZATOP BSS



- Detection of **Influenza virus type A** (H1N1, H3N2) & **type B**
- Result in **8 min.**

Ref. SW40002



BIOSYNEX® RSVTOP BSS



- Detection of **RSV antigens** in nasopharyngeal swabs
- Result in **15 min.**

Ref. 1150006BSS



BIOSYNEX® COVID-Ag BSS. Manufacturer: BIOSYNEX SWISS SA (Switzerland).
Packaging: 25 tests + 25 sterile swabs + 25 pre-filled extraction tubes + 25 dropper caps. Storage: 2 - 30 °C.

BIOSYNEX® INFLUENZATOP BSS. Manufacturer: BIOSYNEX SWISS SA (Switzerland).
Packaging: 20 tests + 20 sterile swabs + 25 pre-filled dropper bottles.
Positive control swabs included for Influenza A and Influenza B
Storage: 2 - 30 °C.

BIOSYNEX® RSVTOP BSS. Manufacturer: BIOSYNEX SWISS SA (Switzerland).
Packaging: 10 tests + 10 sterile swabs + 10 diluent tubes + 2 controls. Storage: 2 - 30 °C.
For further information, please refer to the instructions for use of the product.

A complete range of products for the diagnosis of sexually transmitted infections.

CHEMBIO STAT VIEW® HIV TESTS

WORLD'S SMALLEST SAMPLE VOLUME FOR THE DETECTION OF HIV-1 AND HIV-2 ANTIBODIES IN 15 MINUTES

✓ EASY TO USE

- Quick blood collection; needs only 2.5 µl of blood
- Unique Barrel Design: allows direct transfer of blood from finger to test device

✓ CONVENIENT

- Low hands on time, no bottles, no mixing, no pouring
- Innovative test design prevents accidental spilling

✓ TRUSTED & PROVEN

- Over 13 million tests used professionally
- Sensitivity: 100%; Specificity: 99.8%

✓ FLEXIBLE

Available for self-testing and professional use

STAT VIEW® HIV SELF

Packaging: 1 foil pouch containing 1 test device, 1 buffer cap, 1 lancet, 1 bandage + 1 sterile gauze pad + 1 alcohol swab + 1 test stand. For self-testing.

STAT-VIEW® HIV 1/2

Packaging: 25 pouches containing 1 test device, 1 buffer cap, 1 lancet, 1 bandage + 1 disposable rack. For professional use only.

Manufacturer: CHEMBIO DIAGNOSTICS (USA).



Ref. 60-9508-1

CE 0459

✓ WHO PRE-QUALIFIED

✓ SELF TESTING

✓ PROFESSIONAL USED

TEST FOR HIV IN 3 EASY STEPS WITH STAT-VIEW®

1 Prepare

Fill tip of testing device with 2.5 µl of blood

2 Run

SNAP!
SNAP!
SNAP!

Push hard through the foil until 3 snaps are felt

WAIT
15
MINUTES

3 Read

NEGATIVE

POSITIVE

INVALID

**Refer to the instructions for use provided in the test kit*

DPP® SYPHILIS SCREEN & CONFIRM

SIMULTANEOUS AND SEPARATE DETECTION OF
TREPONEMAL AND NONTREPONEMAL ANTIBODIES FROM
ONLY 10 μ L OF BLOOD, SERUM OR PLASMA

A PARADIGM SHIFT IN CONFIRMATORY TESTING
THAT ENABLES IMMEDIATE TREATMENT

✓ HIGH PERFORMANCE

Excellent performance compared to standard laboratory-based RPR and TP-PA assays*

✓ EASY TO USE

No need of specialized medical training

✓ ACTIONABLE

Provides results in 15 minutes, enabling physicians to make treatment decisions

✓ COST EFFECTIVE

Reduces overtreatment rate, healthcare patient visits, and cost burden

DDP® SYPHILIS SCREEN & CONFIRM ASSAY. Manufacturer: CHEMBIO DIAGNOSTICS (USA).

Packaging: 20 test devices + 20 Microsafe® tubes + 1 vial of running buffer


Storage: 2° to 30° C

For further information, please refer to the instructions for use of the product.

* Data from package insert.



**PROFESSIONAL
USED ONLY**



Ref. 65-9520-0



DPP® HIV-SYPHILIS

COMBINED DETECTION OF HIV 1/2 ANTIBODIES AND
TREPONEMA PALLIDUM ANTIBODIES

✓ IMPROVED PATIENT CARE

Result within 15 minutes

✓ PATIENT-FRIENDLY


Two results from only 10 μ L of whole blood, serum or plasma

✓ RELIABLE

High sensitivity: >99% for HIV
and >94% for *T. pallidum* antibodies*



**PROFESSIONAL
USED ONLY**



Ref. 65-9520-0



DDP® HIV-SYPHILIS. Manufacturer: CHEMBIO DIAGNOSTICS (USA).

Packaging: 20 test devices + 20 Microsafe® tubes + 1 vial of running buffer

Storage: 2° to 30° C

For further information, please refer to the instructions for use of the product.

* Data from package insert.

RAPID TESTS READER

BSX READER®

- ✓ **EASY TO USE**
The touch screen guides you through each step of the result interpretation
- ✓ **FULL TRACEABILITY**
Operator, patient and cassette batch traceability allowing you to meet your quality objectives with peace of mind
- ✓ **STANDARDIZED READING OF RESULTS**
Automated reading and interpretation eliminating user error in the interpretation of results

Product	Reference
BSX READER	5060029
DYMO Labelwriter LW550	502252

WOMEN'S HEALTH

AMNIOQUICK® Card BSX
AMNIOQUICK® DUO+ BSX
PREMAQUICK® BSX

BACTERIOLOGY

TETANOTOP® BSX
EasyTETANOTOP® BSX
Biosynex® S. pneumoniae BSX
Biosynex® L. pneumophila BSX
Biosynex® Strep A BSX

VIROLOGY

INFLUENZATOP® BSX
RSVTOP® BSX
Biosynex® COVID-19 Ag BSX

TOXICOLOGY

DRUGCHECK®+ THC BSX
DRUGCHECK®+10 BSX
DRUGCHECK®+11 BSX

See our product references on the last page.

BSX READER®. Manufacturer: BIOSYNEX SA (France).

URINALYSIS

URITOP® R-50S

Semi-automatic urine test strip reader for qualitative and semi-quantitative detection of 5 to 13 parameters

Compatible with test strips URITOP®+5, URITOP®+7 and URITOP®+11

Product	Reference
URITOP F-50S	FA5640
PRINTER DP-mini	E100

See our product references on the last page.

URITOP® R-50S. Manufacturer: DFLcare (Korea). Distributor: BIOSYNEX SA (France)



Ref. 5060029

CE



Ref. FA5640

CE

CE

BECAUSE TREATMENT ALWAYS STARTS WITH THE CORRECT DIAGNOSIS!

	Reference	Product name	Clinical application	Sample type	CE	Tests/kit
WOMEN'S HEALTH	1090003	AMNIOQUICK®	Detection of IGFBP-1 premature rupture of fetal membranes (strip)	Vaginal swab	1	10
	1090005	AMNIOQUICK® Card	Detection of IGFBP-1 premature rupture of fetal membranes (cassette)	Vaginal swab	1	10
	1090005_BSX	AMNIOQUICK® Card BSX	Detection of IGFBP-1 premature rupture of fetal membranes (for use with BSX READER®)	Vaginal swab	1	10
	1090004	AMNIOQUICK® DUO+	Detection of IGFBP-1 and AFP premature rupture of fetal membranes	Vaginal swab	1	10
	1090004_BSX	AMNIOQUICK® DUO+ BSX	Detection of IGFBP-1 and AFP premature rupture of fetal membranes (for use with BSX)	Vaginal swab	1	10
	1090009	PREMAQUICK®	Prediction of preterm delivery	Vaginal swab	1	10
	1090009_BSX	PREMAQUICK® BSX	Prediction of preterm delivery (for use with BSX READER®)	Vaginal swab	1	10
	SW10003	BIOSYNEX hCG® 20 BSS	Detection of hCG (pregnancy test)	Urine / serum / plasma	1	50
	SW10002	BIOSYNEX hCG® 10 BSS	Detection of hCG (pregnancy test)	Urine / serum / plasma	1	40
BACTERIOLOGY	SW30002	BIOSYNEX® LEGIONELLA PNEUMOPHILA BSS	Detection of L. pneumophila serogroup1 antigen	Urine	1	20
	SW30002_BSX	BIOSYNEX® L.PNEUMOPHILA	Detection of L. pneumophila serogroup1 antigen (for use with BSX READER®)	Urine	1	20
	1030002	BIOSYNEX® S. PNEUMONIAE	Detection of capsule wall polysaccharides CWPS	Urine / CSF	1	20
	1030002_BSX	BIOSYNEX® S. PNEUMONIAE	Detection of capsule wall polysaccharides CWPS (for use with BSX READER®)	Urine / CSF	1	20
	1030005	BIOSYNEX® STREP A	Detection of Strep A antigens (cassette)	Oropharyngeal swab	1	20
	1030005_BSX	BIOSYNEX® STREP A BSX	Detection of Strep A antigens (cassette) - (for use with BSX READER®)	Oropharyngeal swab	1	20
	1030008	BIOSYNEX® C. DIFFICILE GHD	Detection of Clostridium difficile GDH	Feces	2	25
	10025	STREPTATEST®	Detection of Strep A antigens (strip)	Oropharyngeal swab	1	25
	1030011	PYLORITOP® Ag	Detection of Helicobacter pylori antigens	Feces	1	25
	SW30001	PYLORITOP® Ab BSS	Detection of anti-Helicobacter pylori antibodies	WB / serum / plasma	1	25
	1030003	CHLAMYTOP®	Detection of Chlamydia antigens	Endocervical or urethral samples	1	20
	1030016	TETANOTOP®	Detection of antibodies to tetanus antitoxin	WB / serum / plasma	1	20
	65-9525-0	DPP® HIV-Syphilis	Detection of antibodies to HIV-1, HIV-2 and Treponema pallidum	WB / serum / plasma	3	20
	65-9520-0	DPP® Syphilis Screen & Confirm	Detection of treponemal and nontreponemal antibodies	WB / serum / plasma	1	20
VIROLOGY	SW40008	BIOSYNEX® COMBO COVID-FLU BSS	Detection and differentiation of SARS-CoV-2, Influenza A & Influenza B viruses	Nasopharyngeal swab	1	25
	SW40011	BIOSYNEX® COMBO COVID-FLU-RSV BSS	Detection and differentiation of SARS-CoV-2, Influenza Virus (A & B) & RSV viruses	Nasopharyngeal swab	1	25
	SW40006	BIOSYNEX® COVID-19 Ag BSS	Detection of SARS-CoV-2 antigens	Nasopharyngeal swab	1	25
	SW40006_BSX	BIOSYNEX® COVID-19 Ag BSX	Detection of SARS-CoV-2 antigens (for use with BSX READER®)	Nasopharyngeal swab	1	25
	SW40002	BIOSYNEX® INFLUENZATOP	Detection of Influenza A and Influenza B antigens	Nasopharyngeal swab	1	25
	SW40002_BSX	BIOSYNEX® INFLUENZATOP	Detection of Influenza A and Influenza B antigens (for use with BSX READER®)	Nasopharyngeal swab	1	25
	1150006BSS	BIOSYNEX® RSVTOP BSS	Detection of RSV antigens	Nasopharyngeal swab	1	25
	1150006_BSX	BIOSYNEX® RSVTOP BSX	Detection of RSV antigens (for use with BSX READER®)	Nasopharyngeal swab	1	25
	SW40017	BIOSYNEX MONO BSS	Detection of infectious mononucleosis heterophilic antibodies	WB / serum / plasma	1	25
	1150013	IMMUNOQUICK®	Detection of Norovirus, Rotavirus and Adenovirus	Feces	1	20
	1150010	IMMUNOQUICK® NOROVIRUS	Detection of Norovirus	Feces	1	10
	SW40020	BIOSYNEX® COMBO ADENO/ROTA BSS	Detection of Adenovirus and Rotavirus	Feces	1	25
	60-9508-1	STAT-VIEW® HIV Self-test	Detection of HIV-1 and HIV-2 antibodies	WB / serum / plasma	4	1
	60-9520-0	STAT-VIEW® HIV 1/2	Detection of HIV-1 and HIV-2 antibodies	WB / serum / plasma	4	25
	60-9500-0	HIV 1/2 STAT-PAK®	Detection of HIV-1 and HIV-2 antibodies	WB / serum / plasma	3	20
	65-9507-0	DPP® HIV 1/2	Detection of HIV-1 and HIV-2 antibodies	WB / oral fluid / serum / plasma	3	20
	65-9525-0	DPP® HIV-Syphilis	Detection of antibodies to HIV-1, HIV-2 and Treponema pallidum	WB / serum / plasma	3	20
	65-9565-0	DPP® ZCD IgM/IgG	Differential diagnosis of Zika, Chikungunya & Dengue	WB / serum / plasma	1	20
MYCOLOGY	1120001	BIOSYNEX CryptoPS	Semi-quantitative detection of Cryptococcus sp. antigens	WB / serum / plasma / CSF	2	20

BECAUSE TREATMENT ALWAYS STARTS WITH THE CORRECT DIAGNOSIS!

	Reference	Product name	Clinical application	Sample type	CE	Tests/
PARASITOLOGY	5499	PALUTOP® +4 OPTIMA	Detection of 4 Plasmodium species: Plasmodium falciparum (Pf), Plasmodium vivax (Pv), Plasmodium malariae (Pm), Plasmodium ovale (Po)	WB	1	10
	60-9550-0	Chagas STAT PAK®	Detection of anti-Trypanosoma cruzi antibodies	WB / serum / plasma	1	20
	65-9565-0	DPP® ZCD IgM/IgG	Differential diagnosis of Zika, Chikungunya & Dengue	WB / serum / plasma	1	20
BIOCHEMISTRY	1140002_SW	DRUGCHECK®+ AMP	Detection of amphetamines	Urine	1	10
	1140004_SW	DRUGCHECK®+ BUP	Detection of buprenorphine	Urine	1	10
	1140008_SW	DRUGCHECK®+ COC	Detection of cocaine	Urine	1	10
	1140001_SW	DRUGCHECK®+ COT	Detection of cotinine	Urine	1	10
	1140021_SW	DRUGCHECK®+ EDDP	Detection of 2-Ethylidene-1,5-Dimethyl-3,3-Diphenylpyrrolidine	Urine	1	10
	1140022_SW	DRUGCHECK®+ FYL	Detection of fentanyl	Urine	1	10
	1140008_SW	DRUGCHECK®+ MET	Detection of metamphetamines	Urine	1	10
	1140009_SW	DRUGCHECK®+ MOP	Detection of morphine	Urine	1	10
	1140023_SW	DRUGCHECK®+ OXY	Detection of oxycodone	Urine	1	10
	1140011_SW	DRUGCHECK®+ PCP	Detection of phenylcyclidine	Urine	1	10
	1140013_SW	DRUGCHECK®+ THC	Detection of cannabis metabolites	Urine	1	10
	1140013_BSX	DRUGCHECK®+ THC BSX	Detection of cannabis metabolites (for use with BSX READER®)	Urine	1	10
	1140025_SW	DRUGCHECK®+ TML	Detection of tramadol	Urine	1	10
	1140014_SW	DRUGCHECK®+ 5	Detection of AMP, COC, MDMA, MOP, THC	Urine	1	10
	1140015_SW	DRUGCHECK®+ 10	Detection of AMP, BAR, BZO, COC, MDMA, MET, MOP, MTD, TCA, THC	Urine	1	10
	1140015_BSX	DRUGCHECK®+ 10 BSX	Detection of AMP, BAR, BZO, COC, MDMA, MET, MOP, MTD, TCA, THC	Urine	1	10
	1140016_SW	DRUGCHECK®+ 11	Detection of AMP, BAR, BUP, BZO, COC, MDMA, MET, MOP, MTD, TCA, THC	Urine	1	10
	1140016_BSX	DRUGCHECK®+ 11 BSX	Detection of AMP, BAR, BUP, BZO, COC, MDMA, MET, MOP, MTD, TCA, THC	Urine	1	10
	SW50003	BIOSYNEX® D-DIMER BSS	Detection of D-dimers (cut-off 500 ng/mL FEU)	WB	1	10
	SW50001	BIOSYNEX® CRP BSS	Semi-quantitative detection of C-reactive Protein levels	WB / serum / plasma	1	20
URINALYSIS	1110002	HEMOTRUST®	Detection of hemoglobin	Feces	1	25
	1110001	BIOSYNEX® PSA	Detection of Prostate-Specific Antigen (PSA) (cut-off 3 ng/mL)	WB / serum / plasma	1	25
	1110004	BIOSYNEX® Bladder Tumor	Detection of Bladder Tumor Antigen (BTA)	Urine	1	10
	1040007	URITOP+5	Detection of Glucose, Proteins, pH, Blood and Ketones	Urine	1	100
	1040008	URITOP+7	Detection of Glucose, Proteins, pH, Blood, Nitrites, Specific gravity and Leukocytes	Urine	1	100
	1040010	URITOP+11	Detection of Glucose, Proteins, pH, Blood, Ketones, Nitrites, Specific gravity, Leukocytes,	Urine	1	100
	1040011	URITOP+13	Detection of Glucose, Proteins, pH, Blood, Ketones, Nitrites, Specific gravity, Leukocytes,	Urine	1	100

● BSX READER® REFERENCES ● CHEMBIO® PRODUCT

WB = Whole Blood

- (1) CE-IVDD certified
- (2) CE0123-IVDR certified
- (3) CE0459-IVDD certified
- (4) CE0459-IVDR certified on going

In vitro diagnostic medical devices.
Carefully read the instructions for use.

FLY_TDR_EN_V01202310R01

BIOSYNEX

BIOSYNEX S.A.

22 Boulevard Sébastien Brant

67400 ILLKIRCH-GRAFFENSTADEN - France

Phone: +33 3 88 77 57 52 - Fax: +33 3 88 78 76 78

E-mail: customer@biosynex.com

www.biosynex.com